

HETTICH-RACK_L AND HETTICH-RACK_{XL} (LÖWENSTEIN) FOR 16 TUBES *

These flexible racks simplify handling, save time and enable the use of

- tubes of different diameters (15-20 mm)
- tubes of different lengths (100-125 mm or 126-170 mm)




Hettich-Rack_L
for tubes of length 100-125 mm, Handle of gold colour

Cat. No. 60027

- different slant angles** (5° or 20°).

The slant angle of the tubes can be very easily changed without the need for tools.

The racks are of high-quality stainless steel and can be autoclaved.



Hettich-Rack_{XL}
for tubes of length 126-170 mm, Handle of champagne colour

Cat. No. 60028

HETTICH ORBITAL SHAKER (HSM 10) SETS TO FIT IN HETTICH INCUBATORS WITH AND WITHOUT COOLING *

The robust metal casing of the shaker enables it to withstand temperatures up to +65 °C and ensures a long service life.

It is available with different platforms depending on the microbiological application. Each combination of shaker and platform is a complete unit.

TECHNICAL DATA HSM 10

- Shaking frequency: 50 - 250 rpm
- Max. Load: 10 kg
- Dimensions (WxDxH): 410x595x140 mm
- Timer



HSM 10 (**set 1**) including the universal platform (5 tension rollers) for different vessels

Cat. No. 60905



HSM 10 (**set 2**) including platform (anti-slip mat) for a wide range of containers

Cat. No. 60906



HSM 10 (**set 3**) including platform for Erlenmeyer flasks, with a set of 18 clamps of different sizes

Cat. No. 60907

*) More information is available from info@hettichlab.com or can be downloaded from our website www.hettichlab.com

**) The slant of the agar increases the surface area available for growth. It is necessary to use test tubes to hold the agar for micro-organisms that grow very slowly to avoid the cultures drying out.